

Introducing....



MagnaWand™ Cleaning Instruments

A non-magnetic tool designed specifically to clean and disinfect MRI and CAT scan machines. Finally Hospitals will be able to easily clean between patients without losing valuable scanning time.

The **MagnaWand** tool reaches inside the MRI magnet bore or CAT scan gantry to clean and disinfect those difficult to reach places in seconds.

MagnaWand's specifically designed disposable pads can be used with any spray disinfectant or cleaning solution.

Once the area of interest has been cleaned and disinfected, simply eject the disposable pad without touching it with your hands, eliminating the risk of further contamination.

Did You Know...

- More than 50 Million MRI and CAT scan studies performed every year
- The Spread of resistant bacteria in Hospitals HAS TRIPLED SINCE 1980
- MRI and CAT scan machines are NOT routinely cleaned between patients
- Deadly bacteria can survive on hard surfaces for up to an hour



Patent Pending

The MagnaWand weighs 8 ounces and is approx. 36 inches long.



Each MagnaWand cleaning pad weighs 1 ounce and measures 4"x 5" inches.

For More information contact us at:
www.MagnaWand.net

MagnaWand, Inc
505 Lafayette Street
Newtown, PA 18940